Optimierung von Kristallisationsprozessen mittels in-situ Messtechnik

Kolloquium Prozessanalytik 16. November 2007



Dr. Jochen Schöll Dr. Paul Barrett



Overview

Brief introduction of MT PAT tools

- PAT tools in crystallization process development
 - Process understanding
 - Process control
 - Scale-up and production
- Conclusions

Why PAT for Crystallization Processes?

- Crystallization: important unit operation for purification and formulation
- Issues with offline analysis:
 - Sampling & sample preparation
 - Temperature sensitive
 - Sample tends to agglomeration / segregation
 - Provides no process dynamics
 - etc.



Introduction: ReactIR[®] and MonARC[®]

ReactIR – the online in-situ reaction analysis system for the characterization, optimization and control of chemical processes



Introduction: FBRM[®] and PVM[®] Technologies



Lasentec[®] : Award-Winning Innovation



FBRM® Method of Measurement



FBRM® Method of Measurement



PVM® image illustrating the view from the FBRM® Probe Window



Probe detects pulses of Backscattered light

And records measured Chord Lengths

FBRM® Method of Measurement





Thousands of Chord Lengths are measured each second to produce the FBRM® Chord Length Distribution :

The FBRM[®] Chord Length Distribution

Real-time measurement of number and dimension



Dimension (micron)

Overview

Brief introduction of MT PAT tools

PAT tools in crystallization process development Process understanding

- Process control
- Scale-up and production
- Conclusions

Solubility and Metastable Zone Width (MSZW)

- The solubility of a substance in a solvent as a function of temperature and/or solvent composition is the fundamental thermodynamic information that lays the foundation for any crystallization process development.
- The Metastable Zone Width (MSZW) characterizes the ability of a system to remain in a supersaturated state without exhibiting spontaneous nucleation. The MSZW is a kinetic information depending both on nucleation and growth kinetics.
- The fast and reliable determination of solubility and MSZW by combining FBRM and an automatic lab reactor will be shown in the following slides.

The Metastable Zone Width concept



Recommendation for robust crystallization process development: Seed half-way in the MSZ and apply controlled cooling profile

Automated and in parallel

Automated MSZW determination combining FBRM and ALR



How is it done?



Supersaturation controlled crystallization



One step further: crystallization kinetics

Nucleation or growth?

Besides the thermodynamic properties of a system, mainly the **kinetics of the fundamental mechanisms** are determining the properties of a given system. Hereby, nucleation and growth play a significant role, since they are mainly depending on the driving force, i.e. supersaturation



The following example highlights how the combined use of FBRM and PVM allows for the fast and reliable characterization of the underlying fundamental kinetics of a crystallization.

Quantifying crystallization kinetics

- The dependence of nucleation and growth kinetics on the driving force, i.e., supersaturation, can be directly quantified for a given system using in situ monitoring tools, as will be shown in the following example.
- In this example, we will see a seeded drowning-out crystallization:
 - Mother liquor with crystal seeds
 - Linear antisolvent addition at constant T
- To simplify the example, the same experiment was repeated at two different reactor temperatures.



Typical crystallization set-up featuring FBRM and PVM in an automated lab reactor

FBRM trends of Run1

FBRM trends of Run1 at lower temperature (T=25°C)



- Initially small, later steep increase of fine particles (1-5µm)
- Decrease in the number of coarse particles (100-250µm)
- Median is decreasing throughout the process

PVM images of Run1

FBRM trends and PVM images of Run1 at T=25°C



PVM images visualize FBRM measurement data

Distribution comparison of Run1

FBRM CLDs and PVM images of Run1 at T=25°C



- Increase in fines highlighted on the left in No.Wt. CLD
- Slight decrease in size can be seen on the right in Sqr.Wt. CLD
- Quantification of fines increasing using statistics:



FBRM trends of Run2

FBRM trends and PVM images of Run2 at 45°C



- Similar trends as observed in Run 1; however, fines generation less pronounced.
- PVM images confirm the slight increase of fine particles.

Comparison of Run1 and Run2



- Less increase of fine particles at higher temperature.
- Increase instead of decrease in the number of coarse counts at higher temperature.

Comparison of endpoint distribution



- Main results: Run1 features more fines and less coarse material than Run2
- Main effects of higher process temperature:
 - Lower supersaturation at higher temperature generates less nuclei
 - Higher temperature accelerates crystal growth

Both effects lead ultimately to a coarser and more homogeneous product

Quantifying crystallization kinetics

In situ FBRM and PVM are valuable tools for the relative quantification of nucleation and growth kinetics for a series of experiments

- For the determination of the dependence of nucleation and growth kinetics on supersaturation additional information is needed, e.g., the liquid phase concentration as monitored via ATR-FTIR spectroscopy.
- Following publications describe protocols to determine the kinetics of
 - Nucleation: J. Schöll et al., Chem. Eng. Technol., 29(2), 257-264, 2006.
 - Crystal growth: J. Schöll et al., Faraday Discuss., 136, 247-265, 2007
 - **Agglomeration:** C. Lindenberg, J. Schöll, L. Vicum, J. Brozio, and M. Mazzotti, "Precipitation of L-glutamic acid: agglomeration effects", submitted to *Crystal Growth and Design*.
 - Polymorph transformation: J. Schöll et al., *Cryst. Growth Des.*, 6(4), 881-891, 2005.

Once these kinetics are known, they can be used together with Population Balance Modeling to simulate the process under different conditions.

Overview

Brief introduction of MT PAT tools

- PAT tools in crystallization process development
 - Process understanding

Process control

- Scale-up and production
- Conclusions

In-Situ Supersaturation Monitoring & Control of Crystallization Processes

Vince Liotta, Ph.D. & Vijay Sabesan Schering-Plough Research Institute

Real-Time Analytics Users Forum February/March 2005





Experimental Setup Crystallization Monitoring



Supersaturation Control at $S=2=\Delta C$



Temperature & Calculated FTIR Data





Supersaturation Control at $S=1.5 = \Delta C$



Crystal Size Results





Sq Wt Median = $95\mu m$



Sq Wt Median = $107\mu m$

Lower Supersaturation yielded larger Crystals



Overview

- Brief introduction of MT PAT tools
- PAT tools in crystallization process development
 - Process understanding
 - Process control

Scale-up and production

Conclusions

Possible downstream processing issues





PAT from Lab-scale to Production: a Lasentec Approach

André Ridder – 12.06.2007 Group Leader Process Development/PAT Champion Heumann PCS, Pfizer PGM, Feucht, Germany

Why are we Dealing with PAT at Feucht?

Continuous quality monitoring
Supporting health and safety
Saving money



Pictures from the Plant Implementation





Old Process: Crystallization





Old Process: Water Wash



HEUMANN PCS ein Unternehmen der Prizer Gruppe

New One-Pot Procedure





Product Particle Size Distribution



Conclusion and Outlook

- Feucht was capable to optimize the process with Lasentec FBRM[®] on plant scale among the first at Pfizer
- Stable process with an extremely narrow particle size distribution
- Cycle time reduction of 16 hours
- Yield improvement by 10% due to better crystallization
 - Cost reduction by 20%



Overview

- Brief introduction of MT PAT tools
- PAT tools in crystallization process development
 - Process understanding
 - Process control
 - Scale-up and production

Conclusions

Conclusions

- In situ PAT tools such as ATR-FTIR, FBRM and PVM are valuable for the fast and reliable determination of thermodynamic and kinetics properties of a given crystallization process and thus allow for a rigorous process development and optimization.
- PVM allows for the detection of multiple phases, components, shapes, and polymorphs while avoiding any sampling problems. The online microscopic images recorded by PVM can offer unexpected improvements in process understanding.
- FBRM chord length distributions can be used as a fingerprint at the end-point of the particulate process, thus it can help to ensure batchto-batch consistency and simplifies scale-up.
- Finally, FBRM and PVM have an enormous potential for process optimization with respect to downstream issues and product properties and therefore assure fast ROI.

Questions and answers

For further application, product and technology information visit

www.mt.com/crystallization

OR

Contact us at autochem.marketing@mt.com

 OR

+ 1.410.910.8500

Visit <u>www.mt.com/webinars2007</u> for the AutoChem webinar schedule.