



## *Characterising Acrylic Polymers with Thermal Field-Flow Fractionation coupled to Multi Angle Light Scattering*

### **A new application report**

from Postnova Analytics demonstrates how Thermal Field-Flow Fractionation coupled to Multi Angle Light Scattering (TF3-MALS-RI) can be used to characterize different polyacrylate formulations.




### **Acrylic polymers (polyacrylates)**

are widely used in many industrial and consumer applications including paints and coatings. These polymers are often complex combinations of monomer units and the extent of cross-linking or branching determines the polymers' physical properties and therefore determines which application(s) a given formulation will be most applicable for.

### **In Thermal Field-Flow Fractionation (TF3)**

- a separation field is established by applying a temperature gradient perpendicular to the channel. The top channel wall is heated and the bottom wall is cooled, driving polymers towards the cold wall by thermal diffusion. Smaller polymers diffuse into the faster flow profiles of the laminar channel flow and elute to detectors sooner than larger polymers.





### **Characterization of three different polyacrylate formulations**

was undertaken using a Postnova Analytics TF3-MALS-RI system enabling access to their molecular weight distributions as well as their radii of gyration. In addition, plotting the log of molecular weight (log M) versus the log of radius of gyration (log Rg) yields a conformation plot, which can be used to estimate the degree of polymer cross-linking.

### **From the data presented**

in this report, the authors show how TF3-MALS-RI is uniquely able to separate large polyacrylate molecules and provide precise molecular weight and Rg values. By comparison, most chromatography column-based separation techniques would filter out some or all polymer molecules this large, resulting in incorrect determination of their size and molecular weight distributions. Data is also shown that demonstrates how TF3-MALS can be used to elucidate the polymer structure, which can provide insight into why different formulations have different physical properties.

### **For a copy of this new application report**

please visit <https://bit.ly/3h0Pm7U> or contact Postnova Analytics on +49-8191-985-6880 / +44-1885-475007 / +1-801-521-2004 or [info@postnova.com](mailto:info@postnova.com).

**For more information** on the Postnova TF2000 Thermal Field-Flow Fractionation system please visit [www.postnova.com/product/systems/tf2000-thermal-fff.html](http://www.postnova.com/product/systems/tf2000-thermal-fff.html)

### **Postnova Analytics**

Founded in 1997, Postnova Analytics is the inventor and leading international supplier of Field-Flow Fractionation (FFF) systems for markets including biopharmaceuticals, polymers, materials, nanotechnology and environmental sciences. Leveraging its unique and patented modular FFF - Light Scattering Platform, Postnova has been able to provide high performance solutions to a wide range of applications. Today the company's growing product portfolio also includes Flow FFF, Centrifugal FFF, Thermal FFF, Size Exclusion Chromatography (SEC), MALS and DLS. Postnova ([www.postnova.com](http://www.postnova.com)) is located and headquartered in Landsberg am Lech (Germany) and has subsidiaries in the USA, UK and Finland. These offices, in conjunction with a highly qualified and trained distribution network, provide informed local support to customers worldwide.

-----  
**Worldwide HQ**

**Postnova Analytics GmbH**  
Max-Planck-Straße 14  
86899 Landsberg am Lech  
Germany

Tel: +49-8191-985-6880  
Email: [info@postnova.com](mailto:info@postnova.com)  
Web [www.postnova.com](http://www.postnova.com)