



# Neue Methoden der Bildanalyse - Teil 2




Automatisierte  
Partikelform- und -größen Bestimmung mit dem  
FPIA 3000

# 1) FPIA-3000 Produkt Übersicht





- Automatisches Bildanalysesystem für die Messung der Partikelgröße und Partikelform von Suspensionen und Emulsionen
- Größenbereich : 0.8 - 300  $\mu\text{m}$
- Bis zu 300 000 Partikel pro Messung (alle Bilder werden gespeichert)
- Mehr als 20 morphologische Parameter
- “Ein Knopf” Bedienung (SOP Konzept)
- 1 Meßzyklus (Analysenzeit + automatische Spülung) in 2 - 3 min

## 2) Messprinzip, schematischer Aufbau




### 3) Mechanismus der Durchflussmesszelle




- Partikel werden vereinzelt und mit ihrer größten Fläche hin zur CCD Kamera orientiert und liegen im Fokus



# 4) Mechanismus der Durchflussmesszelle



- Gewährleistet, daß alle Partikel im Fokus liegen



# 5) Verwendbare Hüllstrom Flüssigkeiten

- **FPIA-3000** Standard Version
  - Wasser
  - Methanol
  - Ethanol
  - Isopropanol
  - Ethylen Glykol Lösung (25%)
- **FPIA-3000S** lösemittelbeständige Version
  - alles von oben
  - Toluol
  - TCE
  - Azeton
  - Heptan
  - Hexan


## 6) Probenzugabe



Probemenge:  
1ml - 5ml (Standard)








## 7) Bilddatei



# 8) Ergebnisdarstellung 1



File(F) Record(R) Test(T) Settings(S) Maintenance(M) Help(H)



Record ListAnalysis ResultsParticle Image ListFrequency TableMeta DataDetail ResultsMulti ScattergramsGraph(Overlay-View)Graph(Trend-View)

ResetDisplay

SearchSearch

Sort

Filter

	Record Information		Record Information	User Information	Sample Information			Testing Parameters	User-defined Item	Testing Parameters	Analysis Results	
	Record number	Tested	ReAnalyzed	Measurer name	SOP group	Sample name	Sample number	Power Field	Item1	Count method	Density	Pa
14	80	25/11/2004 14:27:03		Service		Magnetic powder	10	LPF->HPF(LPF)		Stop By Time	1088	
15	82	25/11/2004 14:32:33		Service		Magnetic powder	12	LPF->HPF(LPF)		Stop By Time	213	
16	83	25/11/2004 14:34:44		Service		Magnetic powder	13	LPF->HPF(HPF)		Stop By Time	2380	
17	84	25/11/2004 14:38:16		Service		awabreak	14	LPF->HPF(LPF)		Stop By Time	9093	
18	85	25/11/2004 14:40:17		Service		awabreak	15	LPF->HPF(HPF)		Stop By Time	15221	
19	93	25/11/2004 17:01:40		Service	?????	Fiber	2	HPF		Stop By Time	56671	
20	94	25/11/2004 17:10:55		Service	?????	Fiber	3	HPF		Stop By Time	2481	
21	95	25/11/2004 17:23:48		Service	?????	toner	4	HPF		Stop By Time	3033	
22	96	25/11/2004 17:31:56		Service	?????	toner u1min	5	HPF		Stop By Time	3231	
23	97	25/11/2004 17:43:00		Service	?????	abractive	6	HPF		Stop By Time	1720	
24	98	25/11/2004 17:47:35		Service	?????	abractive	7	HPF		Stop By Time	2090	
25	99	25/11/2004 17:53:29		Service	?????	grapfite	8	HPF		Stop By Time	383	
26	100	25/11/2004 17:56:53		Service	?????	grapfite	9	HPF		Stop By Time	5713	
27	101	25/11/2004 18:01:12		Service	?????	grapfite u1min	10	HPF		Stop By Time	13973	
28	102	25/11/2004 18:04:50		Service	?????	grapfite u1min surfactant	11	HPF		Stop By Time	17263	
29	103	25/11/2004 18:13:57		Service	?????	toner	12	HPF		Stop By Time	1985	
30	104	25/11/2004 18:19:55		Service	?????	graphite	13	HPF		Stop By Time	214573	
31	107	25/11/2004 18:44:51		Service	?????	2um	1	HPF		Stop By Time	23212	
32	108	26/11/2004 13:57:41		Service	Latex		1	HPF	0	Stop By Time	13380	
33	110	26/11/2004 16:20:44		Service	Latex	Fiber	1	HPF	0	Stop By Time	1178	
34	112	25/11/2004 17:01:40	27/11/2004 15:15:12	Service	Latex	Fiber	2	HPF	0	Stop By Time	56671	
35	162	29/11/2004 09:59:09		Service	Toner	toner	1	HPF	0	Stop By Time	20500	
36	163	29/11/2004 10:01:15		Service	Toner	toner	2	HPF	0	Stop By Time	20093	
37	164	29/11/2004 10:03:00		Service	Toner	toner	3	HPF	0	Stop By Time	19683	
38	165	29/11/2004 10:04:44		Service	Toner	toner	4	HPF	0	Stop By Time	19694	
39	166	29/11/2004 10:06:29		Service	Toner	toner	5	HPF	0	Stop By Time	19463	
40	170	29/11/2004 13:10:24		Service		algues saines 85	1	HPF	0	Stop By Time	6585	
41	173	29/11/2004 13:26:53		Service		algues senescentes 80	2	HPF	0	Stop By Time	1073	
42	176	29/11/2004 13:40:46		Service		fuseaux2	3	HPF	0	Stop By Time	591	
43	178	29/11/2004 13:50:22		Service		fuseaux 2	4	HPF	0	Stop By Time	207	

# 9) Ergebnisdarstellung 2


Record List Analysis Results Particle Image List Frequency Table Meta Data Detail Results Multi Scattergrams Graph(Overlay-View) Graph(Trend-View)

## Sample Information

Tested	25/11/2004 17:01:40
ReAnalyzed	27/11/2004 15:15:12
ReAnalysis	<input checked="" type="checkbox"/>
Measurer name	Service
ReAnalyzed By	Service
SOP group	Latex
Sample name	Fiber
Sample number	2
Sample type	
Comment	
Power Field	HPF
Count method	Stop By Time
Sheath Liquid	Particle Sheath

Circularity

Max. Distance



Image

Merge

Cancel

Save

## Analysis Parameters

SOP:

Make New

Size :

Lower(%)

Upper(%)

Detailed Se

Shape :

Standard

Extended1

Extended2

Extended3

Extended4

Extended5

Diam:

0.50

um <=

Max. Distance

< 300.00

Shape:



0.200

<=

Circularity

<= 1.000

Cumulative(%)



Standard


User definition

Max. Distance	Mean	13.273
	SD	18.367
	CV	138.38
	Mode	1.328
	Lower%	1.629
Circularity	50 %	6.338
	Upper%	36.346
	Mean	0.653
	SD	0.184
	CV	28.20
Analyzed	Mode	1.000
	Lower%	0.426
	50 %	0.627
	Upper%	0.935
	Density	55274
Analyzed	Small(%)	0.00
	Middle(%)	100.00
	Large(%)	0.00
	Selected(%)	100.00
	Detected(%)	20.459



# 10) Ergebnisdarstellung 3

**Results view:** The results of a single analysis are represented in a 3 graph format – a particle size distribution (green), a particle shape distribution (red) and a scattergram plot of size against shape (blue). The statistical parameters associated with each distribution (mean, mode, lower, median and upper percentile values etc) are also displayed.



**Particle view:** Images of all particles are saved. These images can be viewed and manipulated through the particle viewer. The images can be magnified and sorted on any size or shape parameter allowing the operator to quickly and easily identify anomalies – perhaps agglomerates or the presence of unexpected foreign particles for example.



# 11) Verfügbare morphologische Parameter



## Partikel Durchmesser

Diameter of an equivalent circular area  
Diameter of a circle with equivalent circumference  
Maximum length  
Maximum vertical length  
Major-axis diameter  
Minor-axis diameter  
Vertical Feret's diameter  
Horizontal Feret's diameter  
Vertical diameter of an equivalent area  
Horizontal diameter of an equivalent area  
Maximum vertical diameter  
Maximum horizontal diameter  
Particle circumference  
Edge circumference  
Particle area  
Edge area  
Mean luminosity  
Luminosity distribution value

## Partikel Form


Circularity  
Aspect ratio (width/length)  
Circularity of edge (circumference)  
Circularity of edge (encompassed area)  
Horizontal circumscribing rectangle surface ratio  
Vertical circumscribing rectangle surface ratio  
Main axis circumscribing rectangle surface ratio

## Statistische Parameter

Mean  
Mode  
SD  
CV  
50% value  
Lower (may be set optionally in the 1-49% range)  
Upper (may be set optionally in the 51-100% range)  
Particle number restriction  
Particle rate restriction  
Large particle rate  
Medium particle rate  
Small particle rate

# 12) Zirkularität und Durchmesser

**Zirkularität = Kreisumfang / Partikelumfang**



Fläche des Partikels


Kreis gleicher  
Partikelfläche

Partikelumfang


Durchmesser des Kreise  
gleicher Partikelfläche




# 13) Zirkularität




**1.00**



**0.952**



**0.93**



**0.88**



**0.77**




**0.660**

**(1:5)**


[Wadell's sphericity]

# 14) Messungen an Glaskugeln und Latex




10 µm gemessen:

$10 \pm 0.47$




20 µm gemessen:  $20.43 \pm 1.35$



2.5, 5, 10 & 20 µm

# 15) Pigmente - Vergleich der 10 - 20µm Anteile




Pigment A




Pigment B

# 16) Pigment A



Anteil (Nr.)  
Zirkularität  
< 0,85  
  
= 29%

# 17) Pigment B



Anteil (Nr.)  
Zirkularität  
< 0,85  
= 51%

# 18) FPIA-3000 Anwendungsbereiche

- **Toner** (Rundheit als QC Werkzeug)
- **Keramik**
- **Schmiergelmittel** (z.B. Schleifpapierqualität)
- **Pharma**
- **Biotechnik** (z.B. Proteinkristalle, Wachstumshemmung von Algen )
- **u.a**



# 19) Zusammenfassung

## ① Morphologi G2

- Automatisierte Partikelgrößen - und Form Bestimmung an **Pulvern**
- Messbereich  $0.5\mu\text{m}$  -  $1000\mu\text{m}$
- bis 500.000 Partikel pro Messung
- 10 morphologische Parameter
- 1 Knopf Bedienung

## ② FPIA 3000

- Automatisierte Partikelgrößen - und Form Bestimmung an **Partikeln in Flüssigkeit**
- Messbereich  $0.8\mu\text{m}$  -  $300\mu\text{m}$
- bis 300.000 Partikel pro Messung
- über 20 morphologische Parameter
- 1 Knopf Bedienung

Vielen Dank für Ihre Aufmerksamkeit !