



**Remote sample management unchains scientists from their lab bench.**

**Using high-speed, web-based communications protocols - DP5 Network**

rack decoding software from **Ziath** can be run from a laptop, desktop, phone, or any other network connected device.



**Using DP5 Network**

you can plug all your tube rack readers and scanners into your LAN directly and run them from a PC in your office or somewhere else in the building. This flexibility enables you to have your reader / scanner conveniently near the freezers or compound stores and the results on your desktop or phone.

**DP5 Network**

is fast and secure - it allows you to export tube positions and barcodes in more formats than anyone else - JSON, XLS, CSV, Text, PDF, PNG and Python and send the data in an email to a colleague or your desktop PC.

**While any camera-based reader**

can be upgraded to run DP5 Standard software, only Ziath DataPaq™ Express and DataPaq™ Mirage 2D-barcode rack readers can be currently supplied with LAN adaptors and RJ45 connectors to run the new DP5 Network software.

**For further information**

please visit <https://ziath.com/products/upgrade-to-dp5> or contact Ziath on +44-1223-855021 / +1-858-880-6920 / [info@ziath.com](mailto:info@ziath.com).



**Ziath Ltd.**

specialises in development of innovative instrumentation and software for sample management using 2D Data Matrix bar-coded tubes. These products are widely used in academic, biotech and pharmaceutical organisations. Part of the Consumables and Instruments division of Azenta Life Sciences (Chelmsford, MA), the company is well placed to support customers worldwide.

-----

**Worldwide HQ**

**Ziath Ltd**

Unit 2a, Solopark Trading Estate  
Station Road  
Papisford  
Cambridge CB22 3HB  
UK

Tel: +44 1223 855021

Email: [info@ziath.com](mailto:info@ziath.com)

Web [www.ziath.com](http://www.ziath.com)